

MightyLaser

Laser status

Best performances achievement

Setup at ATF

Experiments at CELIA and results

Upcoming investigations

MightyLaser

Laser status

Best performances achievement


Setup at ATF

Experiments at CELIA and results

Upcoming investigations

Laser amplification system


During january-february 2010, best performances achievements at CELIA with the following setup.



Laser amplification system

During january-february 2010, best performances achievements at CELIA with the following setup.

Output power cs pump power for different signal powers



MightyLaser

Laser status

Best performances achievement

Setup at ATF



Experiments at CELIA and results

Upcoming investigations

Laser amplification system

Amplification efficiency around 34% instead of the 57% measured at the CELIA in Bordeaux, France.

$\lambda=1032 \text{ nm}$
 $\langle P \rangle = 322 \text{ mW}$
 $\langle \tau \rangle = 200 \text{ fs}$
 $f_{\text{rep}} = 178.5 \text{ MHz}$



ATF 04/12/2010

MightyLaser

Laser status

Best performances achievement



Setup at ATF

Experiments at CELIA and results

Upcoming investigations


Laser amplification system

Amplification efficiency around 30% with the 2m connectorised fiber used at ATF.




⇒ A new 2m fiber not connectorised has also been tested with similar efficiency

Residual pump measurement




Transmitted power vs pump power



⇒ The addition of the whole power going out of the fiber doesn't match the pump power even if we consider injection losses.

⇒ This could be explained by photodarkening Inside the fiber

Back to the setup with the best performances achieved



⇒ The efficiency of the main amplifier is still around 30%

MightyLaser

Laser status

Best performances achievement

Setup at ATF

Experiments at CELIA and results

Upcoming investigations

Upcoming investigations

- Be sure that the problem doesn't come from pump injection
=> Use an other pump diode with a 100 μ m fiber instead of 400 μ m in the main amplifier
- Improve the signal injection in the main amplifier
- Measurement of the output power dependence on pump wavelenght
=> Changing the cooling temperature of the pump diode for a given power
- Go to higher output power with a 300W/200 μ m fiber pump diode